

Class 6 Self-Adhesive Indicators

be sure...be ParaSure





Unprocessed indicator Ready for use

Self Adhesive CLASS 6



PASS Items in tray/pack can be safely used

Self Adhesive CLASS 6



PASS Items in tray/pack can be safely used

Successfully processed indicators Colour has changed to teal or darker

Self Adhesive CLASS 6



FAIL Do NOT use items in tray/pack

Self Adhesive CLASS 6



FAIL

Do NOT use items in tray/pack

Unsuccessfully processed indicators Evidence of original purple colour

Self Adhesive CLASS 6



FAIL Do NOT use items in tray/pack

PS1210 3.5 minutes at 134°C PS1220 4.0 minutes at 134°C

PS1230 5.3 minutes at 134°C PS1240 7.0 minutes at 134°C